



VT-1

VISION TESTER

Focus,
To create delicacy.

Each part manifests craftsmanship,
bringing you splendid experience in
different ways.



Complete auxiliary Lens


Built-in auxiliary Lens with comprehensive visual function examination.

First-class Feel

Precise socket design with clear location feedback of every movement.

Multilayer Coating

Multilayer coating, to reduce reflection and improve visual comfort.



Product Features

- Cross optical lens and accurate color filter lens.
- Unique butterfly appearance, dazzling handle wheel, exquisite workmanship and comfortable feel.
- All optical lenses adopt advanced coating process.
- Accurate location when overturning cross cylinders.
- Comprehensive visual function examination and precise and convenient measurement.

Product Parameters

Sphere lens	-19.00D~-+16.75D
	Step0.25D(being 0.12D when 0.12D auxiliary lens is used)
Cylinder lens	0~-6.00D(being 0 ~ - 8.00D when additional lens is used)
	Step 0.25D(being 0.12D when additional lens is used)
Axis	0°~180° ,Step5°
Pupil Distance	50~75mm,Step1mm
Rotary lens	0~20cm/m,Step1cm/m
Prism basal angel	0°~360°
Cross cylinder	±0.25D
Vertex Distance	13.75mm
Size and weight	335(W)×90(D)×320(H)mm/4.5kg

VITOP PTE LTD

Address: 51 Bukit Batok Crescent #08-38, Unity Centre, Singapore 658077

Tel No.: +65 63164749 Email: admin@vitop.com.sg Website: www.vitop.com.sg